

University of Oregon
Zebrafish International Resource Center

TITLE: Embryo Surface Sanitation

Purpose: To minimize the risk of pathogen transmission within a fish population.

Supplies:

- Spray bottle of 70% alcohol (EtOH) or disinfectant solution
- Mesh baskets or strainers (see Fig. 1)
- 5 containers to hold solutions (see Fig. 2)
- Pipettes (for moving eggs/embryos)
- 3-display timer
- Petri dishes
- Squirt bottle containing 0.5X E2 embryo media

Stock Solutions:

0.5X Embryo Media (E2) w/ Methylene Blue:

http://zebrafish.org/zirc/documents/protocols/pdf/Fish_Nursery/E2_solution.pdf

Sodium Hypochlorite (5.65-6%): Fisher Scientific, Catalog # SS290-1

Sodium Thiosulfate Solution: 1 g sodium thiosulfate into 2 L 0.5X E2 embryo medium

Pronase Stock (Concentrate) Solution: 30 mg protease (from *Streptomyces griseus*, Sigma P-6911) into 1 ml 0.5X E2 embryo medium.

- To prepare a 1-2 month supply, add 0.60 grams of protease powder into 20ml nanopure water. Stir VERY gently by hand until completely dissolved. Keep solution on ice throughout preparation and aliquot process.
- For storage, aliquot 0.68 ml of pronase stock solution into 1.5 ml eppendorf tubes, avoiding the creation of bubbles during pipetting, and place into freezer storage at -62°C (-80°F) for long term or -28°C (-20°F) for short term.

Working Solutions:

Load and Rinse Solution (Containers 1 and 5): 0.5X Embryo Media (E2) w/ Methylene Blue

Bleach Solution (Containers 2): 0.59 ml sodium hypochlorite (chemical grade bleach) into 1 L RO water (make fresh daily and do not use expired bleach)

Sodium Thiosulfate Solution (Container 3): Use stock solution listed above

Pronase Solution in Embryo Media (Container 4): 0.68 ml pronase stock solution into 500 ml embryo media (make fresh daily or each time embryos need to be bleached and use within 2 hours of preparation)

Note: All bottles should be changed at least once per week. The pronase bottle should be changed daily.

Procedure:

1. If needed, acclimate your embryos and solutions to approximately 28.5°C prior to bleaching. Ideally, embryos should be bleached between 24 and 30hpf, and should not be bleached if they are beginning to hatch.
2. Sanitize the counter with 70% alcohol (EtOH) or disinfectant.
3. Set up containers with the following solutions (recipes above and see Fig. 2):
 - a. **Container 1 (Load):** Embryo media with methylene blue
 - b. **Containers 2 (Bleach):** Bleach solution
 - c. **Container 3 (Neutralize):** Sodium thiosulfate solution
 - d. **Container 4 (Pronase Chorions):** Pronase solution
 - e. **Container 5 (Rinse):** Embryo media with methylene blue

Starting with container 1, use the smallest volume of solution necessary to cover the eggs; each subsequent container should have a larger volume than the previous container.

Note: The pronase working solution is of low concentration and therefore does not remove the chorions immediately. Allow at least 24 hours for the embryos to hatch following pronase exposure.

4. **Containers 1 and 2 (Load and Bleach Embryos/Eggs):**
 - a. First, remove the unfertilized or dead embryos before bleaching (VERY IMPORTANT).
 - b. Place the required size and number of mesh baskets into container 1 which contains 0.5X embryo media with methylene blue (see Figures). When bleaching more than one stock, label each mesh basket so that you can keep track of the embryos/eggs or stock in each basket.
 - c. Pipette up to 100 un-hatched embryos into the small mesh baskets, 200 into the medium baskets and 300 embryos into the large mesh baskets. **Be sure to change the pipette between stocks.**
 - d. Lift the mesh baskets out of container 1 and place into container 2 (bleach solution).
 - e. Start the timer and **bleach for 10 minutes**. Move the eggs as little as possible or death could occur.
 - f. Transfer baskets to container 3.
5. **Container 3 (Neutralize Bleach in Sodium Thiosulfate Solution):**
 - a. Swirl eggs **for 5 seconds**.
 - b. Transfer baskets to container 4.
6. **Container 4 (Pronase Chorions):**
 - a. Swirl eggs **for 1 minute**.
 - b. Transfer baskets to container 5.
7. **Container 5 (Rinse Embryos):**
 - c. Swirl eggs **for 2 minutes**.

8. Transfer embryos to Petri dishes:
 - a. Rinse the eggs out of the mesh basket into Petri dishes or finger bowls using a squirt bottle of 0.5X E2 embryo media with methylene blue.
 - b. Label the Petri dishes or finger bowls appropriately.




Fig 1. Mesh baskets. Made using hot glue to adhere nylon mesh (~400 microns) to the bottom of PVC pipe. Various sizes can be made to hold larger numbers of eggs.



Fig 2. Embryo sanitation station. 5 containers are set-up for 0.5X embryo media, bleach, sodium thiosulfate, pronase, and rinse solutions. Pre-made large volume solutions are shown in the background.