

Pseudoloma neurophilia

Diagnostic assay

This assay was designed for diagnosis of *Pseudoloma neurophilia* infection of zebrafish (*Danio rerio*) neural tissue. The brain and spinal cord are dissected post-mortem from fresh or frozen tissue and Proteinase K-digested for PCR. Specific primers were designed to amplify a segment of the small subunit rRNA of *P. neurophilia*. Control primers are utilized in each reaction to amplify a segment of zebrafish DNA.

Reagent concentrations in the PCR reaction:

dNTP concentration: 0.3 mM each

P. neurophilia-specific primer concentration (PNA_03 and PNA_04): 0.5 μ M each

zebrafish control primer concentration (POA05 and POA06): 0.125 μ M each

***Pseudoloma neurophilia*-specific primers:**

PNA_03: 5' TGA AAT GTG GTG ACC CGT TTA GG 3'

PNA_04: 5' TCC TTG ACC CAT CCT TCC TGT G 3'

Control primers:

POA05: 5' GCG TCT AGC TTT GCC CTT TGA TG 3'

POA06: 5' CCG TTT TTG AAG ACA TCT GGT CG 3'

PCR program (60_C_40):

1. 94°C for 3 min
2. 94°C for 30 sec
3. **60°C for 40 sec**
4. 72°C for 40 sec
5. Go to step 2 (above) for 39 cycles
6. 72°C for 5 min
7. 8.0°C hold
8. END

Product size: 441 bp for PNA_03 and PNA_04
186 bp for POA05 and POA06